G E R M A N Y

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OPERATING MANUAL

BIOLOGISCHES INVERS-MIKROSKOP
MBL 3200

BIOLOGICAL INVERTED MICROSCOPE
MBL 3200

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
9 Helligkeitsregler

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1 Illumination

2 Insertion for phase contrast equipment

3 Condensor

4 Field diaphragm

5 Photo-/video connection

6 Eyepieces


7 Change between eyepiece head
and photo/video connection

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2. Assembly of the microscope

The microscope and its components, such as the stand, the ocular head, the lenses, the eyepieces and other accessories are separately packed in styrofoam material. Ensure all small parts are removed from the Styrofoam material.

1. Remove the stand from the packaging

2. Head is mounted.

3. Remove the caps from the head.

Set the eyepieces (6) into the eyepiece head.

4. Take off the caps from the lens turret.

Remove the lenses from the plastic boxes and screw them into position in the lens turret.

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5. Insert the power plug into the power source and into the mains supply socket (12).

The basic functions of microscope are ready for use.

3. Operating instructions

1. Turn on the mains switch (8).
2. Put a sample to the stage.
3. Focus with the 10x lens and the coarse adjustment knob onto the sample.
4. Grasp the eyepieces (6), right and left, and move them out or inwards to obtain the correct pupillary distance required for the microscope to provide a normal field of observation. Looking through the right eyepiece with the right eye and focus with the fine adjustment knob (10) to obtain the sharpest image. Then look through the left eyepiece with your left eye and turn the eyesight adjustment to obtain the sharpest image .
5. The depth of field of the image contrast and the resolution capacity of the microscope lens can be adjusted with the a aperture diaphragm (4) built into the condenser (3).
6. Swivel the selected lens into position.
7. Adjust the light intensity with the brightness control (9) on the left side.

Important! Never counter-rotate the right and left adjustment knobs at the same time.

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**4. Exchanging of bulb**

Disconnect the mains plug before every bulb change.
Fingerprints on the glass bulb impair the light efficiency and reduce the life of the bulb.

2. Pull out the faulty bulb and insert the new bulb carefully.

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5. Maintenance

1. Keep the instrument dust-free. Always cover the microscope with the vinyl dust cover. If dust gets onto lenses, blow off with blower or use a fine lens brush. Optical lenses are very soft, and can easily get scratched by dust particles.
2. Never touch the lens. Fingerprints are very difficult to remove. However, if they occur, use soft linen moistened with isopropyl alcohol to remove it.
3. Keep instrument away from moisture. Moisture is the worst enemy of the optical lenses. Once fungus affects the lenses, the instrument cannot be returned to the original condition. Keep in a dry place.
4. Do not drop the microscope. Optical parts are extremely sensitive and may get broken. Although the mechanical parts are made with the highest precision they could be damaged if dropped.
5. In the event of malfunction, your microscope should only be repaired by KRÜSS.

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Power supply: 90-240 VAC

Fuse: 1 A

Illumination: built-in lamp 6V 30W,
brightfield condenser

XY Stage moving: 120 x 80mm
range

Photo-/Video-: Photoadapter
connection videoadapter with eyepiece

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
To convert the microscope to phase contrast usage follow the instructions below:

1. Place the phase contrast lens 20x in the lens turret.
2. Put the phase contrast diaphragm in the insertion (2).
3. Remove one eyepiece from the microscope head.
4. Place the centering telescope in the microscope head at this point.
5. Focus the darker ring pattern of the phase contrast ring (appearing in the field of view) by adjusting the centering telescope eyepiece.
6. Move the lens turret to the highest position by turning the coarse adjustment knobs.
7. Lower the lens turret slowly until a sharp image of the brightly illuminated pattern of the condenser annular stop appears.
8. The brightly illuminated pattern must be adjusted so that it is completely concentric with and overlapped by the darker pattern. Adjust both centring knobs which are on the front of the phase contrast condenser (3). The retardation annulus is intentionally made slightly larger than the image of the annular diaphragm to assure a perfect and complete overlapping and to prevent leakage of light around the retardation annulus. Any such light scattering would affect the quality of contrast.
9. Replace the telescope with the previously removed eyepiece.

The instrument is now ready for use with the phase contrast method.

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9. Photo and Video connections

The Microscope MBL 3200 is equipped with an additional photo tube for connecting a reflex or video camera.

With the aid of the photo adapter, a reflex camera can be connected to the microscope. With the aid of the C-mount adapter , a video camera can be connected to the microscope.