Hinweis:

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ENGLISH

Contents

Foreword	2
Produkt description	3
Accessoires (option)	4
Sicherheitshinweise	4
Safety information	5
Contraindications	6
Recommended reading	7
Before use	8
Changing the laser tip (option)	9
Use when connected to plug power pack (option)	10
Treatment	11
Application	11
Cleaning	11
Technical Data	12
Guarantee	13
Manufacturer and service address	13
Disposal	13

Laser beam

**Do not look into the beam when
using optical instruments.**

Laser class 2 M

CE 1275

We recommend that you read the accompanying instructions for use carefully. It contains important information on what to look for when operating the mini BIOLAS® Softlaser.

Foreword

By buying mini BIOLAS® Softlaser, you have bought a high quality medical laser. The appliance is handy, easy to use and safe, even if you don't have specialist knowledge.

It's completely safe to treat yourself with mini BIOLAS® Softlaser, if you follow the instructions. It's also painless and aseptic and doesn't cause any damaging heating or any side effects.

Product description

The mini BIOLAS® Softlaser is a combination appliance that can be used both for surface treatment and laser acupuncture. For laser acupuncture, the laser beam is concentrated using a lens system so that it comes to a point at the end of the tip.

Further from the tip, the beam diverges strongly and allows the surface of the skin to be treated.


Materials used

Brass-coated or gold-plated.

Please note!

In order to get the best possible use from your mini BIOLAS® Softlaser, do not subject the appliance to extremes of temperature.

Produkt description

Accessories for stationary operation

If you are using the appliance in stationary mode for longer periods of time, we recommend that you use the plug power pack PWR60-150. If the plug power pack does not come with your mini BIOLAS® Softlaser, you can order it separately from our service centre.

Use the mini Biolas® only with original GML plug power pack PWR60-150. Take notice of the instructions for use.

Accessories (option)

you can use a laser tip adapter with an optical wave guide or a massage tip with a lens in place of the laser tip. The adapter or massage tips are available from your BIOLAS trader

Safety informations

Warning! Danger to eyesight.

Never shine the laser beam into the eyes and do not look into the beam, even with optical instruments such as glasses or lenses etc.

The parts of the laser that are related to safety are firmly connected together and with the casing and should not be

separated. If the connections come apart for any reason, the mini BIOLAS® Softlaser should not be operated under any circumstances. Remove the batteries or detach the laser from the plug power pack and send the mini BIOLAS® Softlaser to the manufacturer for repair.

Keep the mini BIOLAS® Softlaser out of the reach of children at all times.

Wounds and bacterial infections may only be treated with surface laser beams at a distance of 0.5 to 1 cm.

Only treat wounds and bacterial infections with a clean laser tip.

Use of nickel

The mounted parts of the mini BIOLAS® Softlaser (head, laser part, laser tip) are hard gold-plated and therefore largely resistant to micro-organisms. For galvanisation, the gold layer has to have a nickel base coat. It is therefore

possible that people with a nickel allergy could suffer an allergic reaction.

Suggestions for use

- Acne / spots
- Arthritis
- Eczema
- Fine lines
- Fever blisters
- Herpes simplex
- Small scars and wounds
- Cold sores
- Muscle tension
- Pustules
- Chapped skin
- Tooth ache
- Strains
- Gout
- 1 and 2 grade burns
- Tiredness
- Inflammation of joints

And much more.

Contraindications:

The metal parts of the mini BIOLAS® Softlaser are hard gold-plated and provide the best possible protection against micro-organisms. For galvanisation, the gold layer has to have a nickel base coat. It is therefore possible that people with a nickel allergy could suffer an allergic reaction.

- Patients with pacemakers: No laser beam around the ribcage or pacemaker
- Pregnancy: No laser beam in the stomach area / lower thorax / certain acupuncture points (!)
- Patients with tumours: No laser beam on the tumour or in the case of suspicion of relapse or metastases
- No laser treatment of people with epiphysis (growth areas of bones) and hollow bones in children
- No laser treatment of the thyroids, ovaries and testicles
- No laser treatment for patients with epilepsy (oligoepilepsy)

Laser treatment of specific points

If you know the position and effects of acupuncture points and medians, you can do laser acupuncture with the mini BIOLAS® Softlaser

Warning!

Avoid treatments in the eye area because of danger to eyesight.

For use in accordance with provisions, please also refer to the instructions for use and safety information.

Recommended reading

You are bound to have benefited from using your mini BIOLAS® Softlaser and want to know more about laser therapy and laser acupuncture. Maybe you are also interested in the scientific theory of the Chinese needle treatment, Chinese teachings on energy and other types of healing.

Here are some examples of relevant books you can find in bookshops:

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GML:

[Advice on gentle health treatments]

You can use these books to learn all about laser therapy: from how to work with lasers and the properties of laser light to the results of investigations, applications, symptoms and diseases.

Wertsch/Schrecke/Küstner:

Akupunkturatlas

[Acupuncture atlas]

WBV, ISBN 3-92-1988-39-X

Gerhard Leibold:

Akupressur

[Akupressure]

Falken-Taschenbuchverlag, ISBN 3-635-6005-9

These publications give you lots of examples of healing, detailed drawings with the location of the acupuncture points that you can obviously be used for treatment with BIOLAS - and they provide clear explanations on energy teaching.

Before use

Inserting and changing the batteries


If the laser beam only appears as a weak light spot, the batteries have run out! Unscrew the head part (on / off switch) of the appliance, remove the used batteries and replace them with new alkaline manganese batteries, type 4LR44, 6 V.

Please ensure that polarity is correct:

+ Pole to the push button switch, - pole to the laser part.

Then screw on the head part again (on / off switch).


Only use leak-proof, branded batteries! "Cheap" batteries may leak and destroy the electrical parts!

Use during battery operation


Remove the cap from the laser part. Place the cap on the head part until it snaps shut. The laser will then switch on automatically.

To switch off, remove the cap from the head part.

Changing the laser tip (option)

If necessary, you can add a laser tip adapter with optical wave guide or massage tip with quartz lens.
For this purpose, unscrew the laser tip and screw the required alternative tip on by hand.

Warning!

The cap can only be replaced if the laser tip has been removed.

Use during plug power pack operation PWR60-150 (option)

Unscrew the head part of the laser.


Remove the batteries.

Screw the adaptor supplied into the laser, in place of the head part.

Stick the plug of the power pack into the socket of the adaptor.

Stick the network appliance in a plug socket with a 230 volt alternating current.

Switch the supply voltage on using the switch.
If the laser beam lights up pale red, the laser is ready for use.

Treatment

You can treat yourself with mini BIOLAS® Softlaser with no medical supervision. However, treatment using mini BIOLAS® Softlaser is not designed to replace essential medical treatment and interventions. So to make the most of your mini BIOLAS® Softlaser, we recommend that you consult your doctor first.

Directions

The minimum duration for the uses described is approx: 3 minutes per treatment.

In acute cases and especially at the beginning of the treatment, the laser beam treatment may be used 5 – 10 times a day.

If symptoms do not lessen and there is no sign of improvement after three to four days of treatment with mini Biolas® Softlase, please seek medical advice.

If medical diagnosis confirms one of the indications given, treatment with mini Biolas® Softlase can be safely continued until the symptoms of the condition have completely disappeared.

No damaging effects are associated with an overdose of laser treatment.

Cleaning

Use a gentle cleaning product for cleaning the coated surface and the metal parts.

Technical data

Operating voltage: 6 volt direct current via an alkaline manganese battery type 4LR44 or plug power pack PWR60-150. Laser diode: Visible red-light laser diode.

Wavelength approx. 650nm

Laser output: < 5 mW

Laser protection class: 2 M

Medical appliance class: II a

Angle of divergence approx. 17 degrees

EC declaration of conformity

In accordance with EC guideline 93/42/EWG appendix VII with V

It was executed a conformity evaluation process, in accordance with appendix V of EC guideline 93/42/EWG and has been proven by LGA InterCert (identification number 1275).

This declaration is valid for all of the abovementioned products bearing the CE identification and that originate from the following factories:

GML German Medical Laser GmbH
Reutlinger Str. 14
75179 Pforzheim

Tel: + 49 (0) 7231 - 13 95 69-0

Fax: + 49 (0) 7231 - 13 95 69-50

Email: info@gm-laser.de

The declaration of Conformity is available on request

Guarantee

The mini BIOLAS® Softlaser is guaranteed for two years, starting from the date of sale (please retain your receipt). Within the guarantee period, we will repair defects free of

charge, if these defects are proven to be attributed to errors in workmanship or materials. Parts that are subject to wear and tear are excluded from the guarantee. In addition, the guarantee claim is not valid for incorrect operation or failure to adhere to the instructions for use, the use of astringent cleaning products, unauthorised interventions or changes to the appliance.

Disposal

Return used batteries to a specialist retailer. This helps recycling and protects the environment.

The appliance should be disposed of appropriately when it reaches the end of its life. Ask your specialist retailer or manufacturer.

Manufacturer and service address

GML German Medical Laser GmbH
Reutlinger Str. 14
75179 Pforzheim

Tel: + 49 (0) 7231 - 13 95 69-0
Fax: + 49 (0) 7231 - 13 95 69-50

Email: info@gm-laser.de

You can find more products from our company on the Internet at www.gm-laser.de