



TESTED AND COMPATIBLE SMARTPHONES AND TABLETS WITH D-HEART DEVICE
(as of 08/10/2025)

iPhone iOS 13 to 26 (not compatible iOS 13.4)

- iPhone 8, 8 plus, X, XS, XR, 11, 11Pro, 12, 12 Pro, 12 Pro Max, 13 mini, 13, 13 Pro, 13 Pro Max, 14, 14 Plus, 14 Pro, 14 Pro Max, 15, 15 Plus, 15 Pro, 15 Pro Max, 16, 16 Plus, 16 Pro, 16 Pro Max, 17, 17 Pro, 17 Pro Max

iPad

- iPad Mini 4, 5th generation
- iPad Air 2
- iPad (5th, 6th, 7th generation)
- iPad Pro (all models)

Android Android 6.0 to 16.0

Smartphone

- Lenovo Moto G5
- Samsung A21s
- Samsung Galaxy A12, A52, A54, A56 5G
- Samsung Galaxy J5, J6, J7, J8
- Samsung Galaxy S5, S6, S7, S8, S9
- Huawei P9, P9 Plus, P10, P10 Plus, P20, P20 Plus
- LG G6, G7, K41s
- Google Nexus 5X, Pixel, Pixel XL, Pixel 2, Pixel 2XL, Pixel 9a
- Honor 8, 90
- Xiaomi Mi 9T, 11 Lite, 12 Lite 5g
- Oppo A53s
- Redmi 9, Note8T, Note 9s, Note 10Pro, Note 12, Note 14Pro

Tablet

- Lenovo tab M10 FHD Plus
- Samsung Galaxy Tab A8
- Samsung S7 FE
- Samsung S7 FE 5G
- Samsung Tab A 2019
- Huawei Mediapad T3 10'
- Mediacom Smartpad 10
- Acer Iconia Tab M10
- Honor Pad 8
- Xiaomi Pad 5



The mobile device, on which the application is run, should have the following characteristics:

iOS:

- Apple smartphone (or tablet) model 6 or later;
- IOS operating system of version 13 or later;
- 1.3 GHz Dual Core processor;
- RAM 1 Gb;
- 8 Mega Pixel rear camera;
- 3264 x 2448 ppi screen resolution;
- Front camera 1.2 mega pixels;
- Video 30 fps.

Android:

- Android smartphone (or tablet) with operating system of version 6 or later;
- 1.2 GHz 4 Core processor; 2 GB RAM;
- 13 Mega Pixel rear camera;
- 4160 x 3120 ppi screen resolution;
- 5 Mega Pixel front camera;
- Video 30 fps.

Please note that if your smartphone is not in this list, the use of D-Heart device may not be supported. If you have any trouble please contact us at support@d-heartcare.com !

Thanks for your collaboration!