


# CE LIGHT mit patentiertem Drehschalter



CE Made in Germany 



**CE 250** mit Glühlampe / with bulb  
**CE 260** mit weißer Power LED / with white power LED

Drehen für Dauerschaltung/  
Turn for steady light

- Robustes Penlight in exklusivem, sehr edlem Design
- Patentierter Drehschalter für Moment- und Dauerschaltung
- Patentierter, sehr stabiler Clip mit spezieller Vertiefung für einen sehr guten Halt am Hemd oder Kittel
- Clip mit Prägung CE und Germany
- Gehäuse aus hochwertigem, lackiertem Aluminium
- Zur Augen- und Rachendiagnostik geeignet
- Robust penlight with exclusive and noble design
- Patented turn-switch for instant and constant connection
- Patented, robust clip with a special recess for a very good hold on the shirt or lab coat
- Clip with stamping CE and Germany
- Housing made of top-quality varnished aluminium
- Suitable for medical diagnostics

Schaltmechanismus / Shifter mechanism	Moment- + Dauerschaltung / Instant connection + Stand-by switch		LED / LED	Glühlampe / Bulb
Maße + Gewicht / Measurements + Weight	L 130 mm Ø 12 mm / 41 g	Leuchtdauer / Burn time	50 h +	5 h +
Material / Material	Aluminium lackiert/Varnished aluminium	Lebensdauer Leuchtmittel / Lifetime illuminant	ca. 100.000 h	Nennwert: 5h
<b>VERPACKUNG/TYPE OF PACKAGING</b>	<b>SCHWARZER SAMTBEUTEL/BLACK VELVET BAG</b>	Lichtstärke / Luminous intensity	22.000 - 25.000 mcd	0,25 A
Batterietyp / Type of battery	AAA, Micro, LR03, R03, AM4	Betriebsspannung / Operating Voltage	3,0 V	2,2 V
Spannungsversorgung / Voltage supply	2 x 1,5 V	<b>INKLUSIVE 2 x ALKALINE-BATTERIEN</b>		

Weitere Etuis / Further boxes: ETU 01; ETU 02; ETU 10; ETU 12

Farben/Colours:

Max. Druckfläche/max. area for printing : L=65 mm, B=30 mm  
 Max. Gravurfläche/ max. area for laser engraving: L=55 mm, B=7mm



silber lackiert  
silver varnished

Ab 500 Stck. lieferbar in jeder Pantone-Farbe / From 500 pcs. available in each Pantone-colour