

# WONDER mit patentiertem Klickschalter



CE Made in Germany ■■■



Drücken für Dauerlicht/  
Push for steady light



Momentschaltung/  
Instant connection

**WO 400** mit Glühlampe / with bulb  
**WO 410** mit weißer Power-LED / with white power-LED

- Patentierter Klickschalter für Moment- und Dauerschaltung
- Exklusives, elegantes Design
- Gehäuse aus hochwertigem, eloxiertem Aluminium
- Ein Glanzlicht für alle Bereiche
- Zur Augen- und Rachendiagnostik geeignet

- Patented click shifter mechanism for instant and constant connection
- Exclusive, classy design
- Top-quality housing made of anodized aluminium
- A highlight for your advertising
- Suitable for medical diagnostic

Schaltmechanismus/Shifter mechanism	Moment- + Dauerschaltung / Instant connection + Stand-by switch		LED / LED	Glühbirne / Bulb
Maße + Gewicht/Measurements + Weight	L 137 mm Ø 12 mm / 36 g	Leuchtdauer/Burn time	50 h +	5 h +
Material/Material	Aluminium eloxiert/Anodized aluminium	Lebensdauer Leuchtmittel/Lifetime illuminant	ca. 100.000 h	Nennwert: 5h
<b>VERPACKUNG/TYPE OF PACKAGING</b>	<b>POLYBEUTEL/POLYBAG</b>	Lichtstärke/Luminous intensity	22.000 - 25.000 mcd	0,25 A
Batterietyp/Type of battery	AAA, Micro, LR03, R03, AM4	Betriebsspannung/Operating Voltage	3,0 V	2,2 V
Spannungsversorgung/Voltage supply	2 x 1,5 V	<b>INKLUSIVE 2 x ALKALINE-BATTERIEN</b>		

Weitere Etuis / Further boxes: ETU 01; ETU 02; ETU 10; ETU 12

## Farben/Colours:

Max. Druckfläche/max. area for printing: L=75 mm B=30 mm  
Max. Gravurfläche/ max. area for laser engraving: L=55 mm B=7 mm



schwarz  
black

anthrazit  
anthracite

silber  
silver

grün  
green

rot  
red

blau  
blue